



# Recombinant Human Group 3 secretory phospholipase A2 (PLA2G3)

<b>Product Code</b>	CSB-EP882150HU-B
<b>Abbreviation</b>	PLA2G3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NZ20
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	G SPALRWYRTS CHLTKAVPGN PLGYLSFLAK DAQGLALIHA RWDahrRLQS CSWEDEPELT AAYGALCAHE TAWGSFIHTP GPELQRALAT LQSQWEACRA LEESPAGARK KRAAGQSGVP GGGHQREKRG WTMPGTLWCG VGDSAGNSSE LGVFQGPDLG CREHDRCPQN ISPLQYNYGI RNYRFHTISH CDCDTRFQQC LQNQHDSISD IVGVAFFNVL EIPCFVLEEQ EACVAWYWWG GCRMYGTVPL ARLQPRTFYN ASWSSRATSP TPSSRSPAPP KPRQKQHLRK GPPHQKGSKR PSKANTTALQ DPMVSPRLDV APTGLQGPQG GLKPQGARWV CRSFRRHLDQ CEHQIGPREI EFQLLNSAQE PLFHCNCTRR LARFLRLHSP PEVTNMLWEL LGTTCFKLAP PLDCVEGKNC SRDPRAIRVS ARHLRRLQQR RHQLQDKGTD ERQPWPSEPL RGPMSFYNNQC LQLTQAARRP DRQQKSWSQ
<b>Source</b>	E.coli
<b>Target Names</b>	PLA2G3
<b>Protein Names</b>	Recommended name: Group 3 secretory phospholipase A2 EC= 3.1.1.4 Alternative name(s): Group III secretory phospholipase A2 Short name= GIII sPLA2 Short name= sPLA2-III Phosphatidylcholine 2-acylhydrolase 3
<b>Expression Region</b>	20-509
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.