



Recombinant Human ADP-ribosylation factor-binding protein GGA3 (GGA3)

Product Code	CSB-EP882151HU-B
Abbreviation	GGA3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZ52
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEAEGESLE SWLNKATNPS NRQEDWEYII GFCDQINKEL EGPQIAVRLI AHKIQSPQEW EALQALTVLE ACMKNCGRRF HNEVGKFRFL NELIKVVSPK YLGDRVSEKV KTKVIELLYS WTMALPEEAK IKDAYHMLKR QGIVQSDPPI PVDRTLIPSP PPRPKNPVFD DEEKSKLLAK LLKSKNPDDL QEANKLIKSM VKEDEARIQK VTKRLHLEEE VNNNVRLLSE MLLHYSQEDS SDGDRELMKE LFDQCENKRR TLFKLASETE DNDNSLGDIL QASDNLSRVI NSYKTIIEGQ VINGEVATLT LPDSEGNQSC SNQGTLDILA ELDTTNSLSS VLAPAPTTPS SGIPILPPPP QASGPPRSRS SSQAEATLGP SSTSNAISWL DEELLCLGLA DPAPNVPPKE SAGNSQWHLL QREQSDLDF SPRPGTAACG ASDAPLLQPS APSSSSSQAP LPPFPAPVV PASVPAPSAG SSLFSTGVAP ALAPKVEPAV PGHHGLALGN SALHHLDALD QLLEEAKVTS GLVKPTTSPL IPTTTPARPL LPFSTGPGSP LFQPLSFQSQ GSPPKGPPELS LASIHVPLES IKPSSALPVT AYDKNGFRIL FHFAKECPPG RPDVLVVVVS MLNTAPLPVK SIVLQAAVPK SMKVKLQPPS GTELSPIFPI QPPAAITQVM LLANPLKEKV RLRKLTFFAL GEQLSTEVGE VDQFPPVEQW GNL
Source	E.coli
Target Names	GGA3
Protein Names	Recommended name: ADP-ribosylation factor-binding protein GGA3 Alternative name(s): Golgi-localized, gamma ear-containing, ARF-binding protein 3
Expression Region	1-723
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the Golgi-localized, gamma adaptin ear-containing, ARF-binding (GGA) family. This family includes ubiquitous coat proteins that regulate the trafficking of proteins between the trans-Golgi network



and the lysosome. These proteins share an amino-terminal VHS domain which mediates sorting of the mannose 6-phosphate receptors at the trans-Golgi network. They also contain a carboxy-terminal region with homology to the ear domain of gamma-adaptins. Alternative splicing of this gene results in two transcript variants.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.