



Recombinant Human Tropomodulin-2 (TMOD2)

Product Code	CSB-YP882161HU
Abbreviation	TMOD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZR1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MALPFQKELE KYKNIDEDEL LGKLSEEEKL QLENVLDDLD PESAMLPAGF RQKDQTQKAA TGPFDREHLL MYLEKEALEQ KDREDFVPFT GEKKGRVFIP KEKPIETRKE EKVTLDPELE EALASASDTE LYDLAAVLGV HLLNNPKFD EETANNKGGK GPVRNVVKGE KVKPVFEEPP NPTNVEISLQ QMKANDPSLQ EVNLNNIKNI PIPTLREFAK ALETNTHVKK FSLAATRSND PVAIAFADML KVNKTLTSLN IESNFITGTG ILALVEALKE NDTLTEIKID NQRQQQLGTAV EMEIAQMLEE NSRILKFGYQ FTKQGPRTTRV AAATKNNDL VRKKRVEADR R
Source	Yeast
Target Names	TMOD2
Protein Names	Recommended name: Tropomodulin-2 Alternative name(s): Neuronal tropomodulin Short name= N-Tmod
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a neuronal-specific member of the tropomodulin family of actin-regulatory proteins. The encoded protein caps the pointed end of actin filaments preventing both elongation and depolymerization. The capping activity of this protein is dependent on its association with tropomyosin. Alternatively spliced transcript variants encoding different isoforms have been described.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.