



Recombinant *Schizosaccharomyces pombe* Uncharacterized RING finger protein C548.05c (SPCC548.05c)

Product Code	CSB-BP882195SXV
Abbreviation	SPCC548.05c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P3U8
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVKRSSHQRQV VLDEDDEENY NNNLDDEKME VLLIPQSNST TFASSDATQM YKKSRISSNS ENKKQIPDTK TLLETFFQKIK KTLECPICTE ALQRPFTTHC GHTYCYECLL NWLKESKSCP TCRQKLYTQP SPAYLVYEIM NVVAASNSGF PLVGINENPA KKQKEVLFDG MFKQEDSHYP RSILVDSGD VLRCARCQWE LENPYHCDHC GFQISDDQDS GREWFWDGEN AESDSSLNGD NTRGGNISTN RAFNNMGHAP ITAVPQDWLR FDEGEFVGS DLESDFSGPG EYDVDDGFID NRATSQLSPV ESDDDFVAPV NGSNGNGITA LDSTDSEEID IMNGFDEERD SGGTNMVSRS ETCYNDGQRY DELRRELADI QNESLDSLNS SSNNSPSHNN IHSRQHPFSS DEDEGNIVTN GTGLRSSQSS SQNRGFDLEP YSEVFTASYP RRQARRTRTI QLDSDEES
Source	Baculovirus
Target Names	SPCC548.05c
Protein Names	Recommended name: Uncharacterized RING finger protein C548.05c
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.