



Recombinant Schizosaccharomyces pombe DASH complex subunit ask1 (ask1)

Product Code	CSB-BP882219SXV
Abbreviation	ask1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P6S5
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNNLEQLERL EQSITLALYE IDANFSKCHR TVTTKILPIV EKYAKNCNTI WDSSKFWKQF FEASANVSLG GVEEPVPVES NPSDQDVMSN STEADLQLHT KNEHLEKRHS FVGKSDFPDA AVQGDNTKNE DFVQSTPKKM DVSLEDISLD DAALTPIPAR MQTPLRKPEN NPHTGRSALL HRVLDTNWQV QVTPREPKNL QSQEVMDIDS SPFVSPSPIS MKMDMPSLND RNSSHLSLF AEFEHESYDS INPSGMSPPK TIQFSPHTMG VGSSQQANER SLSLQRKLET LNDSNDSFVK EEDSWEL
Source	Baculovirus
Target Names	ask1
Protein Names	Recommended name: DASH complex subunit ask1 Alternative name(s): Associated with spindles and kinetochores protein 1 Meiotically up-regulated gene 181 protein Outer kinetochore protein ask1
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.