



Recombinant Schizosaccharomyces pombe Zinc finger protein rsv1 (rsv1)

Product Code	CSB-BP882228SXV
Abbreviation	rsv1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P7D9
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKSYECPFCK RVFHRQEHQV RHIRSHTGEK PFECSPSPCK KRFTRRDELI RHVRTHLRKA LVTPEQTLDV NLHKAPDSKP EGDKSTGQEA DKSQNQSKDG SITDPVQAAV LALSVAYAKP TSVSLSPDLD QAQSKLIEKP RRRSASNATG SLNKKNQDPL RRFSISESAG AAAPTSPSPN SKSPPSNRK NRLQNVLSPI ASNNVPDFNQ NYPTESNPMF LSTPRFQNTN GQRTLTVPVV VWDARQPPTS SRSGPLPSVM YNHFPSVPI PATVNDTSME YYLPNAYPHP TGISLPFYPF DSGIPVSPNI PVSPSSSFVP MYPTTFPSSK PQIVNAPPAP SYFSPGSSFG AQLCANGSTL LRADTKQYGL ALPKITNSNL ISPNQTFNPV KSSVKALPTL EPPSSPSHAT ATSSLHTLFH TAPSRSYD
Source	Baculovirus
Target Names	rsv1
Protein Names	Recommended name: Zinc finger protein rsv1 Alternative name(s): Required for stationary phase viability protein 1
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.