



# Recombinant Schizosaccharomyces pombe Zinc finger protein rsv1 (rsv1)

<b>Product Code</b>	CSB-MP882228SXV
<b>Abbreviation</b>	rsv1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9P7D9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKSYECPFCK RVFHRQEHQV RHIRSHTGEK PFECSPSPCK KRFTRRDELI RHRVTHLRKA LVTPEQTLDV NLHKAPDSKP EGDKSTGQEA DKSQNQSKDG SITDPVQAAV LALSVAYAKP TSVSLSPDLD QAQSKLIEKP RRRSASNATG SLNKKNQDPL RRFSISESAG AAAPTSPSPN SKSPPSNRK NRLQNVLSPI ASNNVPDFNQ NYPTESNPMF LSTPRFQNTN GQRTLTVPVV VWDARQPPTS SRGSLPLSVM YNHFPSVPI PATVNDTSME YYLPNAYPHP TGISLPFYPF DSGIPVSPNI PVSPSSSFVP MYPTTFPSSK PQIVNAPPAP SYFSPGSSFG AQLCANGSTL LRADTKQYGL ALPKITNSNL ISPNQTFNPV KSSVKALPTL EPPSSPSHAT ATSSLHTLFH TAPSRSYD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rsv1
<b>Protein Names</b>	Recommended name: Zinc finger protein rsv1 Alternative name(s): Required for stationary phase viability protein 1
<b>Expression Region</b>	1-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.