



Recombinant Schizosaccharomyces pombe Adenylyl-sulfate kinase (SPAC1782.11)

Product Code	CSB-EP882231SXV
Abbreviation	SPAC1782.11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P7G9
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MATNITFHPG SVTKEERIKF VGHPGMTIWM TGLSASGKST IACALEQYLL QRGVTTYRLD GDNVRFGLNS DLGFSEQDRN ENIRRIGHVA KLFADACVVA VTSFISPYRK DRDQAREFHK KDGLPFIEVY VECPVEVAEQ RDPKGLYKRA RAGEIKEFTG ISAPYEAPIS PEIVVSSHTQ SIEECVEKIV NYLLEKDLIT LK
Source	E.coli
Target Names	met14
Protein Names	Recommended name: Adenylyl-sulfate kinase EC= 2.7.1.25 Alternative name(s): ATP adenosine-5'-phosphosulfate 3'-phosphotransferase Adenosine-5'-phosphosulfate kinase Short name= APS kinase
Expression Region	1-202
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.