



# Recombinant *Schizosaccharomyces pombe* Uncharacterized AAA domain-containing protein C24B10.10c (SPCC24B10.10c)

<b>Product Code</b>	CSB-MP882234SXV
<b>Abbreviation</b>	SPCC24B10.10c
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9P7J5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNPTTKRAIK EIVVYALAFG CSWYAAHKLL STLDPYRQKR QDTVSKSRKR LDEWAGEQVK ELETLELNEY EQIVASQLVL PSEIDVSFDD IGGMDEHVNQ LLQDVLFPK YPEVFDTHGG LLSCPGLLL YGPPGCGKTM LAKALAKQSQ ATFINVSVGL LTKWFGESN KLVDALFTLA RKLEPTIIFI DEIDTFLRQR QRTDHEAMAQ IKAEFMSMWD GLLSGQSRVL VLGATNRPAD IDEAIRRRMP KVFSIPLPNA EQRRKILELY LKKVPLEANF DWNGVVNATA GLSGSYIKEV CRSALSVP RR ELFDKHGNDL EAIKYDIQSG GLRSLKTEDF YHYESLQNV S GIDVE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SPCC24B10.10c
<b>Protein Names</b>	Recommended name: Uncharacterized AAA domain-containing protein C24B10.10c
<b>Expression Region</b>	1-355
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.