



Recombinant Schizosaccharomyces pombe Translin-1 (tsn1)

Product Code	CSB-EP882247SXV
Abbreviation	tsn1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P7V3
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNKSIFIQLQ DQIDKEHSIR EKLTAEVDLL DEKLRVLQLL LANCEQNLEN QEEILEALEI IKSKTRGLAE LASNFPYYKY NGVWDRSIQK VVYLYLLASW TGRLDKSLRP TYLLSLSEV GQILQVPVFP EESTFHLSIE QYLHAVLSLC SELARQSVNS VISGNYHIPF EALNTIQKVH SSFQVLSLKN DSLRRHFDGL KYDLKRSEDV VYDLRIHKLKLV
Source	E.coli
Target Names	tsn1
Protein Names	Recommended name: Translin-1 Alternative name(s): Meiotically up-regulated gene 90 protein Translin homolog
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.