



Recombinant *Cryptococcus neoformans* var. *grubii* Phospholipase B (PLB1), partial

| | |
|--------------------------|---|
| Product Code | CSB-EP882256DYM |
| Abbreviation | PLB1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9P8P2 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Cryptococcus neoformans</i> var. <i>grubii</i> serotype A (strain H99 / ATCC 208821 / CBS 10515 / FGSC 9487) (<i>Filobasidiella neoformans</i> var. <i>grubii</i>) |
| Purity | >85% (SDS-PAGE) |
| Source | E.coli |
| Target Names | PLB1 |
| Protein Names | Recommended name: Phospholipase B EC= 3.1.1.5 Alternative name(s): Lysophospholipase |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |