



Recombinant Mouse Epidermal growth factor-like protein 7 (Egfl7)

Product Code	CSB-EP882551MO-B
Abbreviation	Egfl7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9QXT5
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	EHVYRPSRR VCTVGISGGS ISETFVQRVY QPYLTTCDGH RACSTYRTIY RTAYRRSPGV TPARPRYACC PGWKRTSGLP GACGAAICQP PCGNNGGSCIR PGHCRCVPGW QGDTQCQTDVD ECSTGEASCP QRCVNTVGSY WCQGWEGQSP SADGTRCLSK EGPSVPAPNP TAGVDSMARE EVYRLQARVD VLEQKLQLVL APLHSLASRS TEHGLQDPGS LLAHSFQQLD RIDSLSEQVS FLEEHLGSCS CKKDL
Source	E.coli
Target Names	Egfl7
Protein Names	Recommended name: Epidermal growth factor-like protein 7 Short name= EGF-like protein 7 Alternative name(s): Multiple epidermal growth factor-like domains protein 7 Short name= Multiple EGF-like domains protein 7 NOTCH4-like prot
Expression Region	22-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.