



# Recombinant *Deinococcus radiodurans* Lipoyl synthase (lipA)

<b>Product Code</b>	CSB-MP882734DJA
<b>Abbreviation</b>	lipA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9RWA4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Deinococcus radiodurans</i> (strain ATCC 13939 / DSM 20539 / JCM 16871 / LMG 4051 / NBRC 15346 / NCIMB 9279 / R1 / VKM B-1422)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTQQDPSTKE PKFIKNGIYR KDSVPVREKK PEWLKVTIPT GQVFTEVRKI VKEHRLHTVC EEAMCPNIGE CWSRGATATFM LMGHICTRAC RFCAVDTGNP MGKLDLDEPR SVADSVRLMD LKYVVLTSVD RDDLPDGGAY HFAKTVKAIAK EVNPQTRVEA LTPDFGGNTA CVDLVLDLSDV DTYAQNLETV RRLTHPVRDI RASYDRTLVS LAHAKQARPD VITKTSLMLG LGETREEIRE AMADCRAAGV DVLTFGQYLR PTMHHPVER YISPAEFDEI REEGMQLGFL EVVSGPLVRS SYKAEQIVMD RPGNLPEHLS HLDGGSELTL I
<b>Source</b>	Mammalian cell
<b>Target Names</b>	lipA
<b>Protein Names</b>	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
<b>Expression Region</b>	1-331
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.