



# Recombinant Arabidopsis thaliana Calcium-dependent protein kinase 31 (CPK31)

<b>Product Code</b>	CSB-YP882800DOA
<b>Abbreviation</b>	CPK31
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9S9V0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GCYSSKNLK QSKRTILEKP FVDIGKVYIL GDELGQQQFG ITRKCVKETS GKTYACKTIL KTNLKSREDE EAVKREIRIM KHLSGEPNIV EFKKAYEDRD SVHIVMEYCG GGELFKKIEA LSKDGTKSYSE KEAVEIIRPI VNVVKNCHYM GVMLRDLKPE NFLLSSTDKN ATVKAIDFGC SVFIEEGEVH RKFAGSAYYI APEVLQGKYG KEADIWSAGI ILYILLCGKP PFVTEPEAQM FSEIKSAKID VDSESWKFID VKAKHLVNRM LNRNPKERIS AAEVLGHPWM KDGEASDKPI DGVVLSRLKQ FRDMNKLKQV ALKVIAANLS EEEIKGLKTL FTNIDTDKSG TITLEELKTG LTRLGSNLSK TEVEQLMEAA DVDGNGTIDI DEFISATMHR YRLDRDDHVY QAFQHFDKDN DGHITKEELE MAMKEHGVGD EVSIKQIITE VDTDNDGKIN FEEFRTMMRS GSSLQPQREL LPIK
<b>Source</b>	Yeast
<b>Target Names</b>	CPK31
<b>Protein Names</b>	Recommended name: Calcium-dependent protein kinase 31 EC= 2.7.11.1
<b>Expression Region</b>	2-484
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.