



Recombinant Arabidopsis thaliana Putative F-box protein At3g49520 (At3g49520)

Product Code	CSB-YP882820DOA
Abbreviation	At3g49520
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SCK7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTTISDLPYD LVKEIFSWVP FTSLRAVRST CKTWNALSKN QIFGKKSVAR NQFLELMILD SRVCSLRFDL QKIRNEDEED LVDPSMKQIS IPNNDQVEI SRVYHCDGLL LCIPKDNSSL MLWNPYLGQT KRIRPKNTFH RDDSFALGYN NRNHKILRLN EENESHIDVY DFSSDSWRTV PDDNPYRDEL IYQSGVSLKG NAYFFDREVT TEAEVGTEDL LITGIEDYLR CFDFITTERFG PRLPLPFILP SPSFEYLALS WARDDKLAVL YSHYDTFDII EIWISTKIEP NAVSWSTFLK VDMSLINGLT DDFLIRFEPK SFFIDEEKKV AVLFDTKVTE TCRYQMAYIV GDDGYFKSVN IGVISNSHWI GGELVRSSYV PSELLQLQV
Source	Yeast
Target Names	At3g49520
Protein Names	Recommended name: Putative F-box protein At3g49520
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.