



Recombinant Arabidopsis thaliana Putative FBD-associated F-box protein At3g50710 (At3g50710)

Product Code	CSB-EP882822DOA-B
Abbreviation	At3g50710
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SCQ5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDRISNLSDD LLLKIVSSLP TKDVVVTMLL SKRWKFLWMM VPKLRFDFEF ELEPSYYGRF LKYVDKSMVL NRAQVLETVK FNVGPCCSSE DIATWIRIGM VRNMRELEIS HCEGYFREHR SIKLPKSLYT YEKLEVLKLA STVVLNVPID VCFPSLKSLH LVCVEYKTKK SHRRLLSGCP VLEELVLDKS YNSFHVRSFY VEIPTLQSL S ILDTSGELYG DFTFVVNAPA LKYFNFVDFY GDLCLRDNMP EVVDVNIKVI YRNPCKLLGP LKSVKRLSLC LSPSTTLHNH MEFYQLVHLE LCGDALMWW D LLTWMLQSSP KLQVLKIYEC KCEEHDYLLD PIEEHWEEPS SVPQCLLFHL NIFEWKYYNA GDEEKVVAY ILKNARQLKT ATFSAASYLY PKEERSRELN ELVYMARASS SCQLLLD
Source	E.coli
Target Names	At3g50710
Protein Names	Recommended name: Putative FBD-associated F-box protein At3g50710
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.