



Recombinant Arabidopsis thaliana Mannan endo-1,4-beta-mannosidase 3 (MAN3)

Product Code	CSB-EP882848DOA-B
Abbreviation	MAN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SG94
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	D LGVKSASSDG FVSRKGVQFI LNGKPFYANG FNAYWLAYEA TDSTTRFKIT YVFQNAIHD LTIVRTWGFR DGGYRALQIA PGVYDEKTFQ GLDFAIAEAK RLGIKMIITF VNNYSDFGGR KQYVDWAKNT GQNVSSDDDF YTNPLVKQYY KNHVKTMVNR VNTFTKVEYK DEPTIMGWEL MNEPQCRADP SGKTLTAWMN EMALYVKSVD SKHLLSTGLE GFYGDSSPQR KTSLNPVAAN VLGTDFIANH KLD AIDFASI HSYPDLWFPN LDEKSRLNLL RKWLECHLED AQNILKKPLI LGEFGKPTNT PGYTQAQRDA VFNATFDTIY ESAEKGGPAA GALFWHVISD GMNNFKDPLS IVLSENSTTV NIITEESRKL GLIRGKGKLS KTKI
Source	E.coli
Target Names	MAN3
Protein Names	Recommended name: Mannan endo-1,4-beta-mannosidase 3 EC= 3.2.1.78 Alternative name(s): Beta-mannanase 3 Endo-beta-1,4-mannanase 3 Short name= AtMAN3
Expression Region	20-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.