



Recombinant Arabidopsis thaliana Expansin-B4 (EXPB4)

Product Code	CSB-MP882855DOA
Abbreviation	EXPB4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SHD1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	QNETIDV AGSGTAGVTW YGEPFGAGST GGACGYGSAV ANPPPLYAMVS AGGPSLFNNG KGCCTCYQVV CIGHPACSGS PITVTITDEC PGGPCASEPV HIDLSGKAMG ALAKPGQADQ LRSAGVIRVN YKRAACLYRG TNIVFRMDAG ANPYYSFVV EYENGDDLS NVEIQPAGGS FISMQEMRSA VWKVNSGSAL RGPFNIRLTS GESHKVIVAY NVIPANWKPD ESYRSIVNF
Source	Mammalian cell
Target Names	EXPB4
Protein Names	Recommended name: Expansin-B4 Short name= At-EXPB4 Short name= AtEXPB4 Alternative name(s): Ath-ExpBeta-1.1 Beta-expansin-4
Expression Region	24-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.