



Recombinant Arabidopsis thaliana Putative F-box only protein 9 (FBX9)

Product Code	CSB-BP882862DOA
Abbreviation	FBX9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SI34
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSDLPPDLVE DILSRVPATS LKRLRFTCKQ WNSLFKNRRF TEKHFCKAPK QSHVLLWKDY TVCPMSINLN FSGSSIEFKS VLSLKDSHYN SEQVYIAKVF HCDGLLLCTT KDHRLLVWNP CLGETRWINF ENDYKPYSRF SLGYKNNKSC RSYKILRFWT SYLTPNHIGL RYNIYEFTTD SWRVLIDKVS LNYFLIESEN GVSVFKGNTYW LALDEETNFL LGFDFTMERF KRLCLPSNKN CDTMVLVSVR EEKLSVSHQN FCSSKMDIWM TNRIDSETAM SWKRYFSVEF KILSMCHSLP FCNSFLIDEE KKIVISTVRG NENMVNIIGA YNAYYAEVPV ESTNWPCSPY FVSYPVSLVQ IQQCKSKENN KTSLEM
Source	Baculovirus
Target Names	FBX9
Protein Names	Recommended name: Putative F-box only protein 9
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.