



Recombinant Arabidopsis thaliana Putative F-box/FBD/LRR-repeat protein At2g05300 (At2g05300)

Product Code	CSB-EP882876DOA
Abbreviation	At2g05300
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SJ32
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MCSKGLSQS LKEDRISQLP DPLLTQILNL LPTEEAVKTS VLSTRWRTLW LWVPNLELSF SKFPCFNAFL SFGNLFSDSD RASCVESLKL YLDGNDASYL KPWVDALVKR KIQRLYVRR T LGGYSHVMPL SLYVCDSLVS LRLYRLSLVD AEFVSLPCLK ILRLKDIVFH NETTFFERLVS SCPVLEELKI DVVWWDGNGVY KVHSRSLKRF CFRSSSLRFD ESVPGVVDA PLLCCLIIDD SVSESFVVTD LESNAKFDIS LCFRLWRFDE AKRSIINMFL AWISRVVDMK ICAHTFKLIH RYSESSPLPR FGYMSSLYVT LNASELEWFP IFLRSSPNLK SLILERSGSS HQLSRKAMER VSMSSVPECL LTSLEFVEFK APICGLGPEM MLVWYFLKNS PTLKKLTLPL KSHSTKDDFF KKLLEIPRCS TECEIVIL
Source	E.coli
Target Names	At2g05300
Protein Names	Recommended name: Putative F-box/FBD/LRR-repeat protein At2g05300
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.