



Recombinant Arabidopsis thaliana Putative F-box protein At2g04810 (At2g04810)

Product Code	CSB-BP882877DOA
Abbreviation	At2g04810
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SJ78
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	METCNGLMRK NKRQKVSRSNS DWSKLCPDVL RKIYETLRSP VDSHRAKIVC SNWYSVWKC VKRPLCPLRI IHQGDSPTVG DGNRKLMGFS YNSYCMASG NWLLMVDRCCL KFYIYNLLTK ERIDLPSMES KIRGGQVSFK SKSNNYNFGY LVGPSRKDDI VPDYEAWEV KKSLAVLWVD ETTGDYAVAW TFMRYLFSY IKGDYSWCNL NHNGKSLVLF DMACENNKLY LLTMDHHIKI FNFYGDFTG EQNLYPFNFV EDPSEYVWKR KIVIRRSGEV LIVLSLKKKV QNEEKLLFYI FKMNLESRKW ERVYICGDEM LIFGRGVTAL ALEDLDDGIK SNSIYFVGED VWPDHQEHEH RVSNCGFFDI ATSKIEWPKK IYCFINQNQW FVGGVAY
Source	Baculovirus
Target Names	At2g04810
Protein Names	Recommended name: Putative F-box protein At2g04810
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.