



Recombinant Arabidopsis thaliana DNA repair protein RAD51 homolog 2 (RAD51B)

Product Code	CSB-BP882886DOA
Abbreviation	RAD51B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SK02
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MANKLIGEMG LHTKISNIFA ARNIITAKDA LSMTEFELME LLDVGMKEIR SAISFISEAT SPPCQSARSL LEKKVENEHL SGHLPHTLKG LDDTLCGGIP FGVLTELVGP PGIGKSQFCM KLALSASFPV AYGGLDGRVI YIDVESKFSS RRVIEMGLES FPEVFHLKGM AQEMAGRILV LRPTSLANFT ESIQELKNSI LQNQVKLLVI DSMTALLSGE NKPGAQRQPQ LGWHISFLKS LAEFSRIPIV VTNQVRSQNR DETSQYSFQA KVKDEFKDNT KTYDSHLVAA LGINWAHAVT IRLVLEAKSG QRIIKVAKSP MSPPLAFPFH ITSAGISLLS DNGTELKGGP INTIHARGHS DMINFHGDSC
Source	Baculovirus
Target Names	RAD51B
Protein Names	Recommended name: DNA repair protein RAD51 homolog 2 Short name= AtRAD51B
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.