



Recombinant Arabidopsis thaliana Serine carboxypeptidase-like 7 (SCPL7)

Product Code	CSB-EP882926DOA
Abbreviation	SCPL7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SQX6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	ASIVK SLPGFDGPLP FELETGYIGV GEEEEVQLFY YFIKSERNPQ EDPLLLWLSG GPGCSSISGL LYENGPVNVK IEVYNGTLPS LVSTTYSWTK VSSIIYLDQP VGTGFSYSRT KLVNKPDSG EAKRIHEFLH KWLKGKHQEFS SNPFYVGGDS YCGMVIPALV QEISKGNYVC CKPPINLQGY ILGNPSTENE VDINYRIPYA HGMALISDEL YESMKRICKG KYENVDP RNT KCLKLVGEYQ KCTKRINKAL IITPECVDTS PDCYMYRYLL TTYWANDENV QRALHVNKGS IGEWVRCYFE IPYNHDIKSS VPHYMNNSID GYASLIFSGD HDMEVPYLGT QAWIRSLNYS LIDDWRP WMI GDQIAGYTRT YANKMAFATI KGGGHTPEYK PEESYIMFQR WISGQPL
Source	E.coli
Target Names	SCPL7
Protein Names	Recommended name: Serine carboxypeptidase-like 7 EC= 3.4.16.-
Expression Region	26-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.