



Recombinant Arabidopsis thaliana Ethylene-responsive transcription factor RAP2-12 (RAP2-12)

Product Code	CSB-MP882943DOA
Abbreviation	RAP2-12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SSA8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MCGGAIISDF IPPPRSRRVT SEFIWPDLLK NLKGSKKSSK NRSNFFDFDA EFEADFQGFK DDSSIDCDDD FDVGDVFADV KPFVFTSTPK PAVSAAAECS VFGKKVTGLD GDAEKSANRK RKNQYRGIRQ RPWGKWAAEI RDPREGARIW LGTFKTAEEA ARAYDAAARR IRGSKAKVNF PEENMKANSQ KRSVKANLQK PVAKPNNPS PALVQNSNIS FENMCFMEEK HQVSNNNNNQ FGMTNSVDAG CNGYQYFSSD QGSNSFDCSE FGWSDQAPIT PDISSAVINN NNSALFFEEA NPAKLLKSMD FETPYNTEW DASLDFLNEA AVTTQDNGAN PMDLWSIDEI HSMIGGVF
Source	Mammalian cell
Target Names	RAP2-12
Protein Names	Recommended name: Ethylene-responsive transcription factor RAP2-12 Alternative name(s): Protein RELATED TO APETALA2 12
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.