



Recombinant Arabidopsis thaliana Ethylene-responsive transcription factor RAP2-12 (RAP2-12)

| | |
|--------------------------|--|
| Product Code | CSB-YP882943DOA |
| Abbreviation | RAP2-12 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9SSA8 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MCGGAIISDF IPPPRSRRVT SEFIWPDLLK NLKGSKKSSK NRSNFFDFDA EFEADFQGFK DDSSIDCDDD FDVGDVFADV KPFVFTSTPK PAVSAAAECS VFGKKVTGLD GDAEKSANRK RKNQYRGIRQ RPWGKWAAEI RDPREGARIW LGTFTKAEAA ARAYDAAARR IRGSKAKVNF PEENMKANSQ KRSVKANLQK PVAKPNNPS PALVQNSNIS FENMCFMEEK HQVSNNNNNQ FGMTNSVDAG CNGYQYFSSD QGSNSFDCSE FGWSDQAPIT PDISSAVINN NNSALFFEEA NPAKLLKSMD FETPYNTEW DASLDFLNEA AVTTQDNGAN PMDLWSIDEI HSMIGGVF |
| Source | Yeast |
| Target Names | RAP2-12 |
| Protein Names | Recommended name: Ethylene-responsive transcription factor RAP2-12 Alternative name(s): Protein RELATED TO APETALA2 12 |
| Expression Region | 1-358 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.