



Recombinant Callithrix jacchus Renin (REN)

Product Code	CSB-YP883286CYL
Abbreviation	REN
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9TSZ1
Product Type	Recombinant Protein
Immunogen Species	Callithrix jacchus (White-tufted-ear marmoset)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MALVNITSSV ILTNYMDTQY YGEIGIGTPP QTFKVVFDTG SSNVWVPSSK CSRLYTACVY HKLFDASDSS SYKHNGTELT LRYSTGTVSG FLSQDVITVG GITVTQTFGE VTEMPALPFM LAEFDGVVGM GFSEQAIGKV TPLFDNIISQ GLLKEDVFSF YYNRDSENSQ SLGGQIVLGG SDPQHYEGNF HYINLIRTGL WQIPMKGVSV GSSTLLCEDG CLALVDTGAS YISGSTSSIE KLMEALGAKK RLFDYVVKCN EGPTLPDISF HGGKEYTLT SADYVFQESY SSKKLCTLAI HAMDIPPTG PTWALGATFI RKFYTEFDRG NNRIGFALAR</p>
Source	Yeast
Target Names	REN
Protein Names	Recommended name: Renin EC= 3.4.23.15 Alternative name(s): Angiotensinogenase
Expression Region	61-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>