



Recombinant Bovine Perilipin-2 (PLIN2)

Product Code	CSB-EP883307BO-B
Abbreviation	PLIN2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9TUM6
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	ASVAVEPQL SVVTRVANLP LVSSTYDLVS SAYISRKDQY PYLKSLCEMA EKGMKTITSV AVTSALPIIQ KLEPQIAVAN TYACKGLDRI EEKLPILNQP TNQVVANAKG AMTGAKDAVT TTVTGAKDSV ASTITGVVDR TKGAVTGSVE KTKSVVSGSI NTVLRSRVMQ LMSSGVENAL TKSELLVDQY LPLTKDELEK EAKKVEGFDM VQKPSYYVRL GSLSTKLR SR AYQQALCRVE EAKRKGQETI SQLHSAFNLS ELARKNVHNA NQKIQDAQDK LYLSWLEWKR SIGYDDTDES HCAEHIESRT LAIARNLTQQ LQTMCHTLLS NIQGLPQNIQ DRANHLGVMA GDIYSVFRNA ASFKEVSDGL LASSKGQLQK MKESLDDVMD YLVNNTPLNW LVGPFYPQVT ESESAQAPGT TRRPGRWSRK HPKPVPVSNA EGSQPDDSSS
Source	E.coli
Target Names	PLIN2
Protein Names	Recommended name: Perilipin-2 Alternative name(s): Adipophilin Adipose differentiation-related protein Short name= ADRP
Expression Region	2-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.