



Recombinant Protein SHQ1 homolog (Y48A5A.1)

Product Code	CSB-BP883319CXY
Abbreviation	Y48A5A.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9TYM6
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MLTPVFWITQ DDDALLIRIR APHGNI AELD YDHGDYMFVF TCPPYFLRLH FKQMVEEYGS GNGSVEWKSD EGEFHIKVPK MHKKEHFSNL DMITELLTPS TTHHQPHGNQ LVEEMDDSED DDEGDGSEFL VEQQPAAEPE EPKSDGKIEK FGYGFGWSKF GVIERLRDEI GKIVDILEPE NVEIEKRADK LMEFDWENFD EGRYLADTLE PEEELLAVIS SKFAQKLEIS DEDRTKLKDL KKS KTS AKIN GNDVEIMTSL IDIVFGYCYD QRVNDWESAC ESGWNC AKLS PSL SFFAKFS SVKECLLACT RRALTYPLYR SFQLTQRVIQ DVCHVITAGG GRPALLHILC DLHRIFIESG EFRYILNDLM IADYIFWIQT VPDEILTRI Q TDLTEISGKI DKYDIGWDLE VLEAEAKLAN TQLDS DDEPD N
Source	Baculovirus
Target Names	Y48A5A.1
Protein Names	Recommended name: Protein SHQ1 homolog
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.