



Recombinant Tropomyosin-2 (TPM2)

Product Code	CSB-EP883335PNV
Abbreviation	TPM2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9U5M4
Product Type	Recombinant Protein
Immunogen Species	Podocoryna carnea (Jellyfish)
Purity	≥85% (SDS-PAGE)
Sequence	MSGEEKLGKL RAKLKEITEQ IDDADQKKVE AKHALVDSL A RLEKNEVEVN SAKRRIKLIE KDLEDSSERL KVAEEKLIKV EAEEKKIEEA RNLLAEAESA DDEKMYNIEE EFKESKRTLE SNETKYIEAQ RKGVVISR DV EKTRDKADTL EKRVAVLEQT IASAGESLVE LEEREGESSE REEINEEKLI FLAGQFKESE VRAEAAERSC NVLERNIFET ENEINTWIIQK RKEIEDEMIE MDTVADEPDD E
Source	E.coli
Target Names	TPM2
Protein Names	Recommended name: Tropomyosin-2
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.