



Recombinant Human Set1/Ash2 histone methyltransferase complex subunit ASH2 (ASH2L)

Product Code	CSB-EP883363HU-B
Abbreviation	ASH2L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UBL3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> MAAAGAGPGQ EAGAGPGPGA VANATGAEEG EMKPVAAGAA APPGEGISAA PTVEPSSGEA EGGEANLVDV SGGLETESSN GKDTLEGAGD TSEVMDTQAG SVDEENGRQL GEVELQCGIC TKWFTADTFG IDTSSCLPFM TNYSFHCNVC HHSGNTYFLR KQANLKEMCL SALANLTWQS RTQDEHPKTM FSKDKDIIPF IDKYWECMTT RQRPGKMTWP NNIVKTMSKE RDVFLVKEHP DPGSKDPEED YPKFGLLDQD LSNIGPAYDN QKQSSAVSTS GNLNGGIAAG SSGKGRGAKR KQQDGGTTGT TTKARSDPLF SAQRLPPHGY PLEHPFNKDG YRYILAEPPD HAPDPEKLEL DCWAGKPIPG DLYRACLYER VLLALHDRAP QLKISDDRLT VVGEKGYSMV RASHGVRKGA WYFEITVDEM PPDTAARLWV SQPLGNLQAP LGYDKFSYSW RSKKGTKFHQ SIGKHYSSGY GQGDVLGFYI NLPEDTETAK SLPDITYKDKA LIKFKSYLYF EEKDFVDKAE KSLKQTPHSE IIFYKNGVNQ GVAYKDIFEG VYFPAISLYK SCTVSINFGP CFKYPPKDLT YRPMSDMGWG AVVEHTLADV LYHVETEVDG RRSPPWEP </p>
Source	E.coli
Target Names	ASH2L
Protein Names	Recommended name: Set1/Ash2 histone methyltransferase complex subunit ASH2 Alternative name(s): ASH2-like protein
Expression Region	1-628
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.