



Recombinant Human Mortality factor 4-like protein 1 (MORF4L1)

Product Code	CSB-MP883368HU
Abbreviation	MORF4L1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UBU8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAPKQDPKPK FQEGERVLCF HGPLLYEAKC VKVAIKDKQV KYFIHYSGWN KKS AVRPRRS EKSLKTHEDI VALFPVPEGA PSVHHPLLTS SWDEWVPESR VLKYVDTNLQ KQRELQKANQ EQYAE GKMRG AAPGKKTSG LQQKNVEVKTK KNKQKTPGNG DGGSTSETPQ PPRKKRARVD PTVENEETFM NRVEVKVKIP EELKPWL VDD WDLITRQKQL FYLPAKKNVD SILEDYANYK KSRGNTDNKE YAVNEVVAGI KEYFNVMLGT QLLYKFERPQ YAEILADHPD APMSQVYGAP HLLRLFVRIG AMLAYTPLDE KSLALLLNLYL HDFLKYLAKN SATLFSASDY EVAPPEYHRK AV
Source	Mammalian cell
Target Names	MORF4L1
Protein Names	Recommended name: Mortality factor 4-like protein 1 Alternative name(s): MORF-related gene 15 protein Protein MSL3-1 Transcription factor-like protein MRG15
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.