



Recombinant Human Ubiquitin thioesterase ZRANB1 (ZRANB1)

Product Code	CSB-MP883378HU
Abbreviation	ZRANB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UGI0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSERGIKWAC EYCTYENWPS AIKCTMCRAQ RPSGTIITED PFKSGSSDVG RDWDPSSTEG GSSPLICPDS SARPRVKSSY SMENANKWSC HMCTYLNWPR AIRCTQCLSQ RRTRSPTESP QSSGSGSRPV AFSVDPCEEY NDRNKLNTRT QHWTCSVCTY ENWAKAKRCV VCDHPRPNNI EAIELAETEE ASSIINEQDR ARWRGSCSSG NSQRRSPPAT KRDSEVKMDF QRIELAGAVG SKEELEVDFK KLKQIKNRMK KTDWLFLNAC VGVVEGDAAA IEAYKSSGGD IARQLTAEV RLLNRPSAFD VGYTLVHLAI RFQRQDMLAI LLTEVSQQAA KCIPAMVCPE LTEQIRREIA ASLHQRKGDF ACYFLTDLVT FTLPADIEDL PPTVQEKLFD EVLDRDVQKE LEEESPIINW SLELATRLDS RLYALWNRTA GDCLLDSVLQ ATWGIYDKDS VLRKALHDSL HDCSHWFYTR WKDWESWYSQ SFGLHFSLRE EQWQEDWAFI LSLASQPGAS LEQTHIFVLA HILRRPIIVY GVKYYKSFGR ETLGYTRFQG VYLPLLWEQS FCWKSPIALG YTRGHFSALV AMENDGYGNR GAGANLNTDD DVTITFLPLV DSERKLLHVH FLSAQELGNE EQQEKLLREW LDCCVTEGGV LVAMQKSSRR RNHPLVTQMV EKWLDRYRQI RPCTSLSDGE EDEDEDE
Source	Mammalian cell
Target Names	ZRANB1
Protein Names	Recommended name: Ubiquitin thioesterase ZRANB1 EC= 3.4.19.12 Alternative name(s): TRAF-binding domain-containing protein Short name= hTrabid Zinc finger Ran-binding domain-containing protein 1
Expression Region	1-708
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.