



Recombinant Human LIM domain-containing protein 1 (LIMD1)

Product Code	CSB-EP883381HU-B
Abbreviation	LIMD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UGP4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDKYDDLGLE ASKFIEDLNM YEASKDGLFR VDKGAGNNPE FEETRRVFAT KMAKIHLLQQQ QQQLLQEETL PRGSRGPVNG GGRLGPQARW EVVGSKLTVD GAAKPPLAAS TGAPGAVTTL AAGQPPYPPQ EQRSRPYLHG TRHGSQDCGS RESLATSEMS AFHQPGPCED PSCLTHGDYY DNLSLASPKW GDKPGVSPSI GLSVGSGWPS SPGSDPPLPK PCGDHPLNHR QLSLSSSRSS EGSLGGQNSG IGGRSSEKPT GLWSTASSQR VSPGLPSPNL ENGAPAVGPV QPRTPSVSAP LALSCPRQGG LPRSNSGLGG EVSGVMSKPN VDPQPWFQDG PKSYLSSSAP SSSPAGLDGS QQGAVPGLGP KPGCTDLGTG PKLSPTSLVH PVMSTLPELS CKEGPLGWSS DGSLSVLLD SPSSPRVRLP CQPLVPGPEL RPSAAELKLE ALTQRLEREM DAHPKADYFG ACVKCSKGVF GAGQACQAMG NLYHDTCTFC AACSRKLRGK AFYFVNGKVF CEEDFLYSGF QQSADRCFLC GHLIMDMILQ ALGKSYHPGC FRCVICNECL DGVPFTVDSE NKIYCVRDYH KVLAPKCAAC GLPILPPEGS DETIRVVSMD RDYHVECYHC EDCGLELNDE DGHRCYPLED HLFCHSCHVK RLEKRPSSTA LHQHHF
Source	E.coli
Target Names	LIMD1
Protein Names	Recommended name: LIM domain-containing protein 1
Expression Region	1-676
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.