



Recombinant Human LIM domain-containing protein 1 (LIMD1)

| | |
|--------------------------|--|
| Product Code | CSB-MP883381HU |
| Abbreviation | LIMD1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9UGP4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MDKYDDLGLE ASKFIEDLNM YEASKDGLFR VDKGAGNPE FEETRRVFAT KMAKIHLLQQQ QQQLLQEETL PRGSRGPVNG GGRLGPQARW EVVGSKLTVD GAAKPPLAAS TGAPGAVTTL AAGQPPYPPQ EQRSRPYLHG TRHGSQDCGS RESLATSEMS AFHQPGPCED PSCLTHGDYY DNLSLASPKW GDKPGVSPSI GLSVGSGWPS SPGSDPPLPK PCGDHPLNHR QLSLSSSRSS EGSLGGQNSG IGGRSSEKPT GLWSTASSQR VSPGLPSPNL ENGAPAVGPV QPRTPSVSAP LALSCPRQGG LPRSNSGLGG EVSGVMSKPN VDPQPWFQDG PKSYLSSSAP SSSPAGLDGS QQGAVPGLGP KPGCTDLGTG PKLSPTSLVH PVMSTLPELS CKEGPLGWSS DGSLSVLLD SPSSPRVRLP CQPLVPGPEL RPSAAELKLE ALTQRLEREM DAHPKADYFG ACVKCSKGVF GAGQACQAMG NLYHDTCTFC AACSRKLRGK AFYFVNGKVF CEEDFLYSGF QQSADRCFLC GHLIMDMILQ ALGKSYHPGC FRCVICNECL DGVPFTVDSE NKIYCVRDYH KVLAPKCAAC GLPILPPEGS DETIRVVSMR RDYHVECYHC EDCGLELNDE DGHRCYPLED HLFCHSCHVK RLEKRPSSTA LHQHHF |
| Source | Mammalian cell |
| Target Names | LIMD1 |
| Protein Names | Recommended name: LIM domain-containing protein 1 |
| Expression Region | 1-676 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final |



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.