



Recombinant Human F-box only protein 40 (FBXO40)

Product Code	CSB-MP883385HU
Abbreviation	FBXO40
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UH90
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGKARRSPPG HHRHCEGCFN RHCHIPVEPN TSCLVISCHL LCGATFHMCK EAEHQLLCPL EQVPCLNSEY GCPLSMSRHK LAKHLQVCPA SVVCCSMEWN RWPNVDSSETT LHENIMKETP SEECDTALA LQDQKVLFRS LKMVELFPET REATEEEPTM NGETSVEEMG GAVGGVDIGL VPHGLSATNG EMAELSQEER EVLAKTKEGM DLVKFGQWEN IFSKEHAASA LTNSSASCES KNKNDSEKEQ ISSGHNMVEG EGAPKKKEPQ ENQKQQDVRT AMETTGLAPW QDGVLERLKT AVDAKDYNMY LVHNGRMLIH FGQMPACTPK ERDFVYGKLE AQEVKTVYTF KVPVSYCGKR ARLGDAMLSC KPSEHKAVDT SDLGITVEDL PKSDLIKTTL QCALERELKG HWISESRSID GLFMDFATQT YNFEPEQFSS GTVLADLTAA TPGGLHVELH SECVTRRH NK SSSAFTFTCN KFFRRDEFPL HFKNVHTDIQ SCLNGWFQHR CPLAYLGCTF VQNHFRPPGQ KAKVIYSQEL KTFAIKPEVA PELSEGRKNN HLLGHGGKSQ NSLTSLPLEI LKYIAGFLDS VSLAQLSQVS VLMRNICATL LQERGMVLLQ WKKKRYSHGG TSWRVHREIW QFSSLFSKIK SWEFNEVTSM SEHLKSCPFN IVEHKTDPI L TSMCQPREQ ARESLVSTFR IRPRGRYVS
Source	Mammalian cell
Target Names	FBXO40
Protein Names	Recommended name: F-box only protein 40 Alternative name(s): Muscle disease-related protein
Expression Region	1-709
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.