



Recombinant Human LIM domain and actin-binding protein 1 (LIMA1)

Product Code	CSB-BP883388HU
Abbreviation	LIMA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHB6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MESSPFNRRQ WTSLSLRVTA KELSLVNKNK SSAIVEIFSK YQKAAEETNM EKKRSNTENL SQHFRKGTLT VLKKKWENPG LGAESHTDSL RNSSTEIRHR ADHPPAEVTS HAASGAKADQ EEQIHPR SRL RSPPEALVQG RYPHIKDGED LKDHOSTESKK MENCLGESRH EVEKSEISEN TDASGKIEKY NVPLNRLKMM FEKGPTQTK ILRAQSR SAS GRKISENSYS LDDLEIGPGQ LSSSTFDSEK NESRRNLELP RLSETSIKDR MAKYQAAVSK QSSSTNYTNE LKASGG EIKI HKMEQKENVP PGPEVCITHQ EGEKISANEN SLAVRSTPAE DDSRDSQVKS EVQQPVHPKP LSPDSRASSL SESSPPKAMK KFQAPARETC VECQKTVYPM ERLLANQQVF HISCFRCSYC NNKLSLGT YA SLHGRIYCKP HFNQLFKSKG NYDEGFGHRP HKDLWASKNE NEEILERPAQ LANARETPHS PGVEDAPIAK VGVLAASMEA KASSQKEKED KPAETKKLRI AWPPPT ELGS SGSALEEGIK MSKPKWPPED EISKPEVPED VDLDLKKLRR SSSLKERSRP FTVAASFQST SVKSPKTVSP PIRKGWSMSE QSEESVGG RV AERKQVENAK ASKKNNGNVGK TTWQNKESKG ETGKRSKEGH SLEMENENLV ENGADSD EDD NSFLKQQSPQ EPKSLNWSSF VDNTFAEEFT TQNQKSQDVE LWEGEVVKEL SVEEQIKRNR YYDEDEDEE
Source	Baculovirus
Target Names	LIMA1
Protein Names	Recommended name: LIM domain and actin-binding protein 1 Alternative name(s): Epithelial protein lost in neoplasm
Expression Region	1-759
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.