



# Recombinant Human Septin-9 (41526)

<b>Product Code</b>	CSB-EP883392HU-B
<b>Abbreviation</b>	41526
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UHD8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKKSYSGGTR TSSGRLRRLG DSSGPALKRS FEVEEVETPN STPPRRVQTP LLRATVASST QKFQDLGVKN SEPSARHVDS LSQRSPKASL RRVELSGPKA AEPVSRRETEL SIDISSKQVE NAGAIGPSRF GLKRAEVLGH KTPEPAPRRT EITIVKPQES AHRRMEPPAS KVPEVPTAPA TDAAPKRVEI QMPKPAEAPT APSPAQTLEN SEPAPVSQLQ SRLEPKPQPP VAEATPRSQE ATEAAPSCVG DMADTPRDAG LKQAPASRNE KAPVDFGYVG IDSILEQMRR KAMKQGFEFN IMVVGQSGLG KSTLINTLFK SKISRKSVQP TSEERIPKTI EIKSITHDIE EKGVRMCLTV IDTPGFGDHI NNENCWQPIM KFINDDQYKEY LQEEVNINRK KRIPDTRVHC CLYFIPATGH SLRPLDIEFM KRLSKVVNIV PVIKADTLT LEERVHFKQR ITADLLSNGI DVYPQKEFDE DSEDRLVNEK FREMIPFAVV GSDHEYQVNG KRILGRKTKW GTIEVENTH CEFAYLRDLL IRTHMQNIKD ITSSIHFEAY RVKRLNEGSS AMANGMEEKE PEAPEM
<b>Source</b>	E.coli
<b>Target Names</b>	SEPT9
<b>Protein Names</b>	Recommended name: Septin-9 Alternative name(s): MLL septin-like fusion protein MSF-A Short name= MLL septin-like fusion protein Ovarian/Breast septin Short name= Ov/Br septin Septin D1
<b>Expression Region</b>	1-586
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.