



Recombinant Human Ubiquilin-2 (UBQLN2)

Product Code	CSB-BP883393HU
Abbreviation	UBQLN2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHD9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AENGESSGP PRPSRGPAAA QGSAAAPAEP KIIKVTVKTP KEKEEFVAVPE NSSVQQFKEA ISKRFKSQTD QLVLIFAGKI LKDQDTLIQH GIHDGLTVHL VIKSQNRPQG QSTQPSNAAG TTTTSASTPR SNSTPISTNS NPFGLGSLGG LAGLSSLGLS STNFSELQSQ MQQQLMASPE MMIQIMENPF VQSMLSNPDL MRQLIMANPQ MQQLIQRNPE ISHLLNNPDI MRQTLEIARN PAMMQEMMRN QDLALS NLES IPGGYNALRR MYTDIQE PML NAAQE QFGGN PFASV GSSSS SGEGTQPSRT ENRDPLPNPW APPPATQSSA TTSTTTSTGS GSGNSSSNAT GNTVAAANYV ASIFSTPGMQ SLLQQITENP QLIQNMLSAP YMRSMMQSLS QNPDLAAQMM LNSPLFTANP QLQEQMRPQL PAFLQQMQNP DTLSAMSNPR AMQALMQIQQ GLQTLATEAP GLIPSFTPGV GVGVLGTAIG PVGPVTPIGP IGPIVPFTPI GPIGPIGPTG PAAPPGSTGS GGPTGPTVSS AAPSETTSPT SESGPNQQFI QQM VQALAGA NAPQLPNPEV RFQQQLEQLN AMGFLNREAN LQALIATGGD INAAIERLLG SQPS
Source	Baculovirus
Target Names	UBQLN2
Protein Names	Recommended name: Ubiquilin-2 Alternative name(s): Chap1 DSK2 homolog Protein linking IAP with cytoskeleton 2 Short name= PLIC-2 Short name= hPLIC-2 Ubiquitin-like product Chap1/Dsk2
Expression Region	2-624
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes an ubiquitin-like protein (ubiquilin) that shares high degree of similarity with related products in yeast, rat and frog. Ubiquilins contain a N-terminal ubiquitin-like domain and a C-terminal ubiquitin-associated domain. They physically associate with both proteasomes and ubiquitin ligases; and thus, are thought to functionally link the ubiquitination machinery to the proteasome to affect in vivo protein degradation. This ubiquilin has also been



shown to bind the ATPase domain of the Hsp70-like Stch protein.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.