



Recombinant Human SERTA domain-containing protein 1 (SERTAD1)

Product Code	CSB-MP883406HU
Abbreviation	SERTAD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHV2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLSKGLKRKR EEEEEKEPLA VDSWWLDPGH TAVAQAPPAV ASSSLFDLSV LKLHHSLLQQS EPDLRHLVLV VNTLRRRIQAS MAPAAALPPV PSPPAAPSVA DNLLASSDAA LSASMALLE DLSHIEGLSQ APQPLADEGP PGRSIGGAAP SLGALDLLGP ATGCLLDDGL EGLFEDIDTS MYDNELWAPA SEGLKPGPED GPGKEEAPEL DEAELDYLMD VLVGTQALER PPGPGR
Source	Mammalian cell
Target Names	SERTAD1
Protein Names	Recommended name: SERTA domain-containing protein 1 Alternative name(s): CDK4-binding protein p34SEI1 Short name= SEI-1 Transcriptional regulator interacting with the PHD-bromodomain 1 Short name= TRIP-Br1
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.