



Recombinant Human Nuclear receptor-binding protein (NRBP1)

Product Code	CSB-EP883410HU-B
Abbreviation	NRBP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHY1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SEGESQTVL SSGSDPKVES SSSAPGLTSV SPPVTSTTSA ASPEEEEESE DESEILEESP CGRWQKRREE VNQRNVPID SAYLAMDTTE GVEVVWNEVQ FSERKNYKLQ EEKVRVAFDN LIQLEHLNIV KFHKYWADIK ENKARVIFIT EYMSSGSLKQ FLKKTCKNHK TMNEKAWKRW CTQILSALSY LHSCDPPIIH GNLTCDTIFI QHNGLIKIGS VAPDTINNHV KTCREEQKNL HFFAPEYGEV TNVTTAVDIY SFGMCALEMA VLEIQNGGES SYVPQEAISS AIQLLEDPLQ REFIQKCLQS EPARRPTARE LLFHPALFEV PSLKLLAAHC IVGHQHMIPE NALEEITKNM DTSAVLAEIP AGPGREPVQT LYSQSPAEL DKFLEDVRNG IYPLTAFGLP RPQQPQQEEV TSPVVPPSVK TPTPEPAEVE TRKVVLMOQN IESVEEGVKH HLTLLKLED KLNRLHSCDL MPNENIPELA AELVQLGFIS EADQSRLTSL LEETLNKFNF ARNSTLNSAA VTVSS
Source	E.coli
Target Names	NRBP1
Protein Names	Recommended name: Nuclear receptor-binding protein
Expression Region	2-535
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.