



Recombinant Human Protein fat-free homolog (FFR)

Product Code	CSB-BP883424HU
Abbreviation	FFR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UID3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>AAAAAAGPS PGSGPGDSPE GPEGEAPERR RKAHGMLKLY YGLSEGEAAG RPAGPDPLDP TDLNGAHFDP EVYLDKLRRE CPLAQLMDSE TDMVRQIRAL DSDMQTLVYE NYNKFISATD TIRKMKNDFR KMEDEMDRLA TNMAVITDFS ARISATLQDR HERITKLAGV HALLRKLQFL FELPSRLTKC VELGAYGQAV RYQGAAQAVL QQYQHLP SFR AIQDDCQVIT ARLAQQLRQR FREGGSGAPE QAECVELLLA LGEPAEELCE EFLAHARGRL EKELRNLEAE LGPSPAPDV LEFTDHGGSG FVGGLCQVAA AYQELFAAQG PAGAEKLA AF ARQLGSR YFA LVERRLAQEQ GGGDNSLLVR ALDRFHRLR APGALLAAAG LADAATEIVE RVARERLGHH LQGLRAAFLG CLTDVRQALA APRVAGKEGP GLAELLANVA SSILSHIKAS LAAVHLFTAK EVSFSNKPYF RGEFCSQGV R EGLIVGFVHS MCQTAQSFCD SPGEKGGATP PALLLLLSRL CLDYETATIS YILTLTDEQF LVQDQFPVTP VSTLCAEARE TARRLLTHYV KVQGLVISQM LRKSVETRDW LSTLEPRNVR AVMKRVVEDT TAIDVQVGLL YEEGVRKAQS SDSSKRTFSV YSSSRQQGRY APSYTPSAPM DTNLLSNIQK LFSERIDVFS PVEFNKVSVL TGIKISLKT LLECVRLRTF GRFGLQQVQV DCHFLQLYLW RFVADEELVH LLLDEVVASA ALRCPDPVPM EPSVVEVICE RG</p>
Source	Baculovirus
Target Names	VPS51
Protein Names	Recommended name: Protein fat-free homolog Alternative name(s): Another new gene 2 protein
Expression Region	2-782
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.