



Recombinant Human Dual specificity protein phosphatase 13 (DUSP13)

Product Code	CSB-MP883436HU
Abbreviation	DUSP13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UII6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSLQKQDLR RPKIHGAVQA SPYQPPTLAS LQRLLWVRQA ATLNHIDEVW PSLFLGDAYA ARDKSKLIQL GITHVVNAAA GKFQVDTGAK FYRGMSLEY Y GIEADNPFF DLSVYFLPVA RYIRAALSVP QGRVLVHCAM GVSRSATLVL AFLMICENMT LVEAIQTVQA HRNICPNSGF LRQLQVLDNR LGRETGRF
Source	Mammalian cell
Target Names	DUSP13
Protein Names	Recommended name: Dual specificity protein phosphatase 13 EC= 3.1.3.16 EC= 3.1.3.48 Alternative name(s): Dual specificity phosphatase SKRP4 Testis- and skeletal-muscle-specific DSP
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Members of the protein-tyrosine phosphatase superfamily cooperate with protein kinases to regulate cell proliferation and differentiation. This superfamily is separated into two families based on the substrate that is dephosphorylated. One family, the dual specificity phosphatases (DSPs) acts on both phosphotyrosine and phosphoserine/threonine residues. This gene encodes different but related DSP proteins through the use of non-overlapping open reading frames, alternate splicing, and presumed different transcription promoters. Expression of the distinct proteins from this gene has been found to be tissue specific and the proteins may be involved in postnatal development of specific tissues. A protein encoded by the upstream ORF was found in skeletal muscle, whereas the encoded protein from the downstream ORF was found only in testis. In mouse, a similar pattern of expression was found. Multiple alternatively spliced transcript variants were described, but the full-length



sequence of only some were determined.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.