



Recombinant Human RWD domain-containing protein 2A (RWDD2A)

Product Code	CSB-MP883450HU
Abbreviation	RWDD2A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UIY3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSASVKESLQ LQLLEMMLF SMFPNQGEVK LEDVNALTNI KRYLEGTREA LPPKIEFVIT LQIEEPKVKI DLQVTMPHSY PYVALQLFGR SSELDRHQQL LLNKGLTSYI GTFDPGELCV CAAIQWLQDN SASYFLNRKL VYEPSTQAKP VKNTFLRMWI YSHHIYQQDL RKKILDVGKR LDVTGFCMTG KPGIICVEGF KEHCEEFWHT IRYPNWKHIS CKHAESVETE GNGEDLRLFH SFEELLLEAH GDYGLRNDYH MNLGQFLEFL KKHKSEHV FQ ILFGIESKSS DS
Source	Mammalian cell
Target Names	RWDD2A
Protein Names	Recommended name: RWD domain-containing protein 2A
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.