



Recombinant Human Protein FAM83C (FAM83C)

Product Code	CSB-MP883580HU
Abbreviation	FAM83C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BQN1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFGGPGPGVL GAQGMAGPLR GRVEELKLPW WRESSPLVLR HSEAARLAAD ALLERGEAAY LRVISEEREL PFLSALDVDY MTSHVRGGPE LSEAQQQEAS GPDRLSLLSE VTSGTYFPMA SDIDPPDL DL GWPEVPQATG FSPTQAVVHF QRDKAKNIKD LLRFLFSQAH TVVAVVMDIF TDMELLCDLM EASSRRGVPV YLLLAQEHLR HFLEMCYKMD LNGEHLPNMR VRSTCGDTYC SKAGRRFTGQ ALEKFVLIDC EQVVAGSYSF TWLCSQAHTS MVLQLRGRIV EDFDREFRCL YAESQPVEGF CGGEDPLSPR ALRPPPVALA FRPDVPSPTS SLPSSTSLSS IKQSPLMGRS SYLALPGGGD CSDTGVVSSS LGPARREASG QPSLHRQLSD PNHGSPPLY RANLGKLGAY PWSQSSPALN HNSTSPLTLA VGSPLLPRSR PLLQFHRGAP ALSRFPENGL PGSQEPSPLR GRWVPGTTLE TVEEKEKKAS PSQSRGQLDL LVPFPRAREV GDPDSGVTPN SGPLRPGEQA PEDRRLSPSQ ADSQLDLLSR ALGTGGAPEL GSLRPGDRAL EDRRLSLNQS RGQSDLLMQY PKAQGSRVPL ETNSSARPAR RAPDERRQTL GHSQDLITK FGPFRGEGPG PNGLPISPA RTAGAGSGDE KRLTLGHSKL DLITKYHQLH GARQGTEPGG PKGGHLNGGN SDLVRDEKRL TLGHSKLDLI TKYNKSKFKQ LRSRFES
Source	Mammalian cell
Target Names	FAM83C
Protein Names	Recommended name: Protein FAM83C
Expression Region	1-747
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.