



Recombinant Human Haloacid dehalogenase-like hydrolase domain-containing protein 3 (HDHD3)

Product Code	CSB-MP883592HU
Abbreviation	HDHD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSH5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAHRLQIRLL TWDVKDTLLR LRHPLGEAYA TKARAHGLEV EPSALEQGFR QAYRAQSHSF PNYGLSHGLT SRQWWLDVVL QTFHLAGVQD AQAVAPIAEQ LYKDFSHPCT WQVLDGAEDT LRECRTRGLR LAVISNFDRR LEGILGGLGL REHFDFVLTS EAAGWPKPDP RIFQEALRLA HMEPVVAAHV GDNLYLCDYQG PRAVGMHSFL VVG PQALDPV VRDSVPKEHI LPSLAHLLPA LDCLEGSTPG L
Source	Mammalian cell
Target Names	HDHD3
Protein Names	Recommended name: Haloacid dehalogenase-like hydrolase domain-containing protein 3
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.