



Recombinant Human Protein fuzzy homolog (FUZ)

Product Code	CSB-MP883593HU
Abbreviation	FUZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BT04
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGEEGTGGTV HLLCLAASSG VPLFCRSSRG GAPARQQLPF SVIGSLNGVH MFGQNLEVQL SSARTENTTV VWKSFHDSIT LIVLSSEVGI SELRLERLLQ MVFGAMVLLV GLEELTNIRN VERLKKDLRA SYCLIDSFLG DSELIGDLTQ CVDCVIPPEG SLLQEALSGF AEAAGTTFVS LVVSGRVVAA TEGWWRLGTP EAVLLPWLVG SLPPQTARDY PVYLP HGSPT VPHRLLTLTL LPSLELCLLC GPSPPLSQLY PQLLERWWQP LLDPLRACLP LGPRALPSGF PLHTDILGLL LLHLELKRCL FTVEPLGDKE PSPEQRRRL RNFYTLVTST HFPPEPGPPE KTEDEVYQAAQ LPRACYLVLG TEEP GTG VRL VALQLGLRRL LLLLSPQSPT HGLRSLATHT LHALTPLL
Source	Mammalian cell
Target Names	FUZ
Protein Names	Recommended name: Protein fuzzy homolog
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.