



Recombinant Human DNA-directed RNA polymerase III subunit RPC7-like (POLR3GL)

Product Code	CSB-MP883594HU
Abbreviation	POLR3GL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BT43
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASRGGRGR GRGQLTFNVE AVGIGKGDAL PPPTLQPSPL FPPLEFRPVP LPSGEEGEYV LALKQELRGA MRQLPYFIRP AVPKRDVERY SDKYQMSGPI DNAIDWNPDW RRLPRELKIR VRKLQKERIT ILLPKRPPKT TEDKEETIQK LETLEKKEEE VTSEEDKEE EEEEEKEEEE EYDEEEHEE ETDYIMSYFD NGEDFGGDSD DNMDEAIY
Source	Mammalian cell
Target Names	POLR3GL
Protein Names	Recommended name: DNA-directed RNA polymerase III subunit RPC7-like Short name= RNA polymerase III subunit C7-like Alternative name(s): DNA-directed RNA polymerase III subunit G-like
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.