



Recombinant Human Acidic leucine-rich nuclear phosphoprotein 32 family member E (ANP32E)

Product Code	CSB-MP883598HU
Abbreviation	ANP32E
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BTT0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEMKKKINLE LRNRSPEEVT ELVLDNCLCV NGEIEGLNDT FKELEFLSMA NVELSSLARL PSLNKLRLKLE LSDNIISGGL EVLAEKCPNL TYLNLSGSKI KDLSTVEALQ NLKNLKSLLD FNCEITNLED YRESIFELLQ QITYLDGFDQ EDNEAPDSEE EDDDEDGEDD EEEEEENEAGP PEGYEEEEEEE EEEEEDEDE DEDEAGSELG EGEEVGLSY LMKEEIQDEE DDDDYVEEGE EEEEEEGGL RGEKRKRDAE DDGEEEDD
Source	Mammalian cell
Target Names	ANP32E
Protein Names	Recommended name: Acidic leucine-rich nuclear phosphoprotein 32 family member E Alternative name(s): LANP-like protein Short name= LANP-L
Expression Region	1-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.