



# Recombinant Human Telomerase Cajal body protein 1 (WRAP53)

<b>Product Code</b>	CSB-MP883605HU
<b>Abbreviation</b>	WRAP53
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BUR4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKTLETQPLA PDCCPSDQDP APAHPSPHAS PMNKNADSEL MPPPPERGDP PRLSPDPVAG SAVSQELREG DPVSLSTPLE TFGSPSELS PRIEEQELSE NTSLPAEEAN GSLSEEEANG PELGSGKAME DTSGEPAED EGD TAWNYSF SQLPRFLSGS WSEFSTQPEN FLKGCKWAPD GSCILTNSAD NILRIYNLPP ELYHEGEQVE YAEMVPVLRM VEGDTIYDYC WYSLMSSAQP DTSYVASSSR ENPIHIWDAF TGELRAS FRA YNHLDELTA A HSLCFSPDGS QLFCGFNRTV RVFSTARPGR DCEVRATFAK KQGQSGIISC IAFSPAQPLY ACGSYGRSLG LYAWDDGSPL ALLGGHQGGI THLCFHPDGN RFFSGARKDA ELLCWDLRQS GYPLWLSL GRE VTTNQRIYFD LDPTGQFLVS GSTSGAVSVW DTDGPGNDGK PEPVLSFLPQ KDCTNGVSLH PSLPLLATAS GQRVFPEPTE SGDEGEELGL PLLSTRHVHL ECRLQLWWCG GAPDSSIPDD HQGEKGQGGT EGGVGELI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	WRAP53
<b>Protein Names</b>	Recommended name: Telomerase Cajal body protein 1 Alternative name(s): WD repeat-containing protein 79 WD40 repeat-containing protein encoding RNA antisense to p53
<b>Expression Region</b>	1-548
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.