



Recombinant Human Uncharacterized protein C20orf195 (C20orf195)

Product Code	CSB-MP883611HU
Abbreviation	C20orf195
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BVV2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSTHVAGLGL DKMKLGNPQS FLDQEEADDQ QLLEPEAWKT YTERRNALRE FLTSDLSPHL LKRHHARMQL LRKCSYYIEV LPKHLALGDQ NPLVLPSALF QLIDPWKFQR MKKVGTAQTK IQLLLLGDLL EQLDHGRAEL DALLRSPDPR PFLADWALVE RRLADVSAVM DSFLTMMVPG RLHVKHRLVS DVSATKIPHI WLMLSTKMPV VFDRKASAAH QDWARLRWFV TIQPATSEQY ELRFRLDPR TQQECAQCGV IPVAACFDV RNLLPNRSYK FTIKRAETST LVYEPWRDSL TLHTKPEPLE GPALSHSV
Source	Mammalian cell
Target Names	FNDC11
Protein Names	Recommended name: Uncharacterized protein C20orf195
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.